

## Safe Use of Water Condensers and Their Alternatives in Research Labs

---

Water cooling is essential in many chemical manipulations and processes, particularly in solvent distillation and purifications, refluxing, and prolonged heating experiments. However, the effectiveness of cooling relies on proper setup and secure connections. Unsafe (loose) tubing connections and water leaks can lead to costly damage and dangerous conditions in the lab. Laboratory personnel can minimize these risks and ensure smooth and safe chemical manipulations by taking necessary safety measures and providing hands-on training on the proper setup, monitoring, and maintenance of cooling systems.

Alternatives to traditional water condenser setups, such as Findensers,<sup>TM</sup> water pumps and recirculating chillers, offer considerable advantages in water conservation, safety, and operational efficiency. These systems are designed to reduce reliance on a continuous water supply, which helps mitigate risks like flooding and water-related accidents in laboratory environments.

On our campus, like most research institutions, there have been several instances of water leaks from condensers. When a cooling system fails, it can cascade into a series of unexpected consequences, from chemical contamination to hazardous reactions to widespread building floods. To prevent similar incidents, it's critical to establish rigorous safety protocols, ensure routine inspections, and train personnel on the proper handling of both condensers and chemical manipulations. By incorporating preventive measures, regular checks, and emergency preparedness into daily laboratory operations, you can reduce the likelihood of such incidents occurring again and mitigate the impact of any unforeseen issues that arise.

### Connecting Tubing to Glass Condensers

Softening the tubing end is an important step to prevent both glass cuts and damage to the condenser. By softening the Tygon tubing's end with a heat gun before inserting it into the condenser port, it is easier to achieve a tight, secure fit while maintaining the integrity of the tubing and condenser. (**Fig. 1**). Never insert tubing with too much force; doing so will damage the condenser and can risk breaking glass and puncturing the skin. It is recommended that you use cut-resistant gloves when doing this or any other task involving force and glassware. Use doubled-back copper wire (wrap-around, twist over, wrap back, twist to tightness), a well-secured plastic zip tie (pull tight with pliers), or a clamp to hold the connection once the tubing has attached to the condenser port; otherwise, it may pop off once the water pressure within the building rises. Due to decreased water use in the building at night, water pressure typically rises. It is important to choose a ridged connector port rather than a smooth one (**Fig. 1**).




Fig. 1: Safe insertion and removal of a Tygon tubing into a glass connector

## Safe Operation of Water Condensers


- Provide comprehensive hands-on training for all lab personnel on properly handling water condensers used for chemical manipulations and emergency procedures. Ensure they are familiar with the potential risks of specific systems and chemicals.
- Establish clear SOPs for the use and maintenance of condensers, including step-by-step instructions on setting up the cooling water system, checking connections, and safely handling the equipment. Regular reviews and updates to these procedures can help prevent accidents.
- Over time, even glass or metal condensers can suffer from wear, corrosion, or cracks. Regular visual inspections can help detect any early signs of damage that could compromise the system's performance and unwanted events in the lab.



Fig. 2. A copper and stainless steel tubing clamp




- Copper wiring, a plastic zip tie or a stainless-steel clamp (**Fig. 2**) should be installed on connections to prevent the tubes from falling apart or bursting the condenser in case of a sudden increase in water pressure. Be cautious when using a clamp, as over-tightening the clamp may crack the condenser port.
- Copper wire is good for hose barbs if you **a)** use a suitable gauge wire and **b)** don't twist so much that you slip over a hump of the barb as you tighten it. We recommend using 14-gauge copper wire for this purpose, which is likely thick enough to provide a strong, secure hold without being too stiff to manipulate. Twisting the wire too much can create hitches when tightening it, as it could cause the wire to slip over the barb or create bumps that make securing it harder. Bending the twisted portion of the wire over after cutting it short is a good safety measure to reduce the risk of injury, ensuring there are no sharp ends that could cause cuts or snags.
- Use tubing with an appropriate diameter and thickness made of Tygon or rubber on all water-cooled systems. Over time, Tygon can become more brittle or less pliable, which can impact its performance in applications such as water transport, or laboratory settings. So, researchers should replace Tygon water tubing every few years.
- Some glass condensers feature easy-to-remove hose barbs, which are handy when setting up and removing distillations routinely (**Fig. 3**). After the barb has been wired or clamped once, any screw-equipped condenser can be readily connected.



**Fig. 3.** A condenser with easy-to-remove hose barbs

- Ensure the water enters the lower end of the condenser and exits from the upper end. This maximizes the heat exchange efficiency by allowing the water to flow in the same direction to the vapor, promoting better cooling and prevent air gaps.
- Avoid using joints, like two-way plastic connectors, between the condenser and the chilled water source (**Fig. 4**). These connectors are prone to detachment, which can disrupt the

flow and function of the system. Instead of connectors, it's recommended to use a single length of Tygon tubing without any joints. Position the reaction assembly as close to the water source as possible. This minimizes the need for long tubing runs and helps to avoid the use of connectors altogether.




**Fig. 4.** A proper reaction set-up inside a chemical fume hood showing correct clamping of Tygon tubing

- Make sure the chemical fume hood cup sink is not clogged with unwanted materials every time you use a water condenser.
- Do not rely on Parafilm® for securing Tygon tubing connections.
- For unattended reactions, laboratory lights should be left on, and “*Safe Operation Card*” signs should be posted near the experiment identifying the nature of the experiment and hazardous substances in use.
- Maintain a modest flow of chilled water throughout the reaction and check the condition of the chilled water tap's valve regularly. Because the facility is not utilized as frequently, the water pressure is low during the day and increases at night.
- Running cold water through the condenser slowly also reduces water use.
- When using a hot plate or an oil bath, make sure that the tubing is not touching the hot surface of the plate. The tubing will melt, and water will flood in the lab.
- Ensure all condensers, tubing, and cooling systems are regularly inspected for wear, leaks, or damage. Set up a scheduled checklist for daily, weekly, and monthly checks, depending on the equipment's usage. Tubing is very affordable; replace it at the first signs of wear.

## Alternatives to Water Condensers in Research Labs

There are a few substitutes for water condensers, such as Findensers™, localized water pumps and recirculating chillers (**Fig. 5 & 6**).




**Fig. 5.** A Findenser

Findensers™ are specialized condensers designed to efficiently reduce the amount of coolant used for chemical reactions, especially when larger-size flasks or high-boiling solvents are involved. Their design uses no external fluid and improves cooling efficiency by combining glass and aluminum components. The Findenser's main component is a standard glass condenser tube. This glass tube provides the vapor with the surface area to condense as it cools. A metal jacket with fins surrounds the glass tube. These fins are crucial because they expand the surface area of the aluminum jacket, improving its ability to disperse heat. The sealed space between the internal glass condenser and the exterior aluminum jacket contains a minimal amount of fluid (**Table 1**). Due to this design, Findensers are particularly effective when dealing with large volumes of high-boiling solvents.

However, Findensers™ may not be suitable for low-boiling solvents like diethyl ether (boiling point ~34-36 °C) and pentane (boiling point ~35-36 °C). If the condenser cannot maintain a sufficiently low temperature, the vapor may escape without effectively condensing, leading to solvent loss. Since the fluid is sealed within the system and doesn't need to be constantly replenished, water consumption is eliminated, which can be especially important in settings where water conservation is a priority and to prevent water flooding from water-cooled condensers.

**Localized water pumps** allow for the use of small volumes of water contained within a closed system. For example, a recirculating pump can be placed in a bucket of ice water to circulate cool water through the condenser. These pumps are typically submersible due to their internal sealing. These pumps significantly reduce the risk of flooding in the laboratory (**Fig. 6**). Traditional water condensers often rely on open water flow, which, if improperly connected or uncontrolled, can lead to overflows or leaks. In contrast, the localized system is typically sealed, and any excess coolant can be contained or recycled within the loop.



**Fig. 6.** A small water pump

The localized water pump system can succeed in short duration experiments requiring precise temperature control because of its small-scale, self-contained design. This is accomplished by controlling the temperature of the contained water. This is particularly true for procedures that need precise temperature control, such as product purifications and solvent distillations. By directly cooling parts, the pump helps maximize performance without overcooling or wasting energy.

**Recirculating chillers** avoid nearly all the problems of water condensers, Findensers™, and simple water pumps but require considerable initial investment (thousands of dollars). In these systems a working fluid (usually water or water/ethylene glycol mixtures) is circulated through tubing to a glassware condenser. The electric chiller unit can be set to a fixed temperature, and it is often possible to set it below ambient (down to -20 °C is a common class of chiller). In the event of a leak in the system, risk for flooding is minimal due to lower total volume, but risks for the chemistry remain (only Findensers™ are truly passive). In an academic setting, appropriately sizing the chiller cooling capacity rarely poses a problem. Take care that operating below ambient

with regular tubing can lead to ice accumulation on the tubing; aiming for a temperature around 0 to  $-10\text{ }^{\circ}\text{C}$  is usually optimal.

**Table 1: Comparison Table Findenser™ vs Water Condenser**

<b>Feature</b>	<b>Findenser™</b>	<b>Water condenser</b>
<b>Efficiency</b>	Less efficient (lower heat capacity of air).	More efficient (higher heat capacity of water).
<b>Size</b>	Findensers™ are surprisingly heavy, and this difference needs to be considered in setting up an experiment. It cannot be dried easily for anhydrous reactions because of the size and fins.	More compact for similar performance. Standard technique in organic syntheses.
<b>Water Consumption</b>	None.	Requires constant water supply from source (possibility of water flooding).
<b>Applications</b>	Small-scale or air-cooled systems	Small, large-scale industrial or lab systems.
<b>Maintenance</b>	Low; Findensers™ cannot be cleaned using an acid wash because they contain metal.	Higher (water treatment and scale management).
<b>Environmental Impact</b>	None.	High consumption of water.
<b>Failure Mode</b>	Medium.	High.

## Training

---

The laboratory personnel should be trained to use any cooling system before conducting any work with chemical manipulations. The principal investigator (PI) or lab manager must provide training to his or her laboratory personnel specific to proper reaction set-up, the hazards involved with chemical manipulation, work area decontamination, and emergency procedures.

Personnel should be trained on how to correctly set up and maintain cooling systems, including water condensers or their alternatives. This includes ensuring secure connections, checking for leaks, and confirming proper water flow rates or the operation of alternative cooling methods (e.g., Findensers™ or water pumps).

**NOTE:** EH&S is available to assist in the education and training of personnel concerning safe work practices.

## References and Important Links

1. Lehman, J. W. Multiscale Operational Organic Chemistry: A Problem-Solving Approach to the Laboratory Course; Prentice-Hall: Upper Saddle River, NJ, **2002**; pp 585 –587.
2. Pappenhagen, L. A. An Improved Reflux Condenser. *J. Chem. Educ.* **1929**, 6, 1618.
3. Solomon, S.; Brook, B.; Rutkowsky, S.; Bennet, J. Using Ice Cooled Condensers in Chemistry Laboratory. *J. Chem. Educ.* **2003**, 80, 299 –301.
4. Waterless Condensers for the Teaching Laboratory: An Adaptation of Traditional Glassware. Baum, E. W.; Esteb, J. J.; and Wilson, A. M. *J. Chem. Educ.* **2014**, 91, 1087–1088.
5. University of Toronto, Department of Chemistry, Flood prevention:  
<https://www.chemistry.utoronto.ca/health-and-safety-guide/flood-prevention>
6. Findensers™: <https://www.radleys.com/range/findenser-super-air-condenser/>
7. Water pumps: <https://www.amazon.com/Aquarium-Pumps-Filters-boxtech>.
8. Recirculating chillers: <https://julabo.us/filter-groups/recirculating-chillers/>.
9. Personnel experience; John Berry, <https://chem.wisc.edu/staff/berry-john-f/>;  
<https://berry.chem.wisc.edu/sops/>.
10. Personnel experience, Prof. Dan Weix, <https://weixgroup.chem.wisc.edu/>.
11. Personnel experience; Rob Toreki, [Safety Emporium](https://www.safetyemporium.com/).
12. Personnel experience; Prof. Craig Merlic, <https://www.chemistry.ucla.edu/directory/merlic-craig/>.
13. VWR Cut Resistant Gloves: [https://us.vwr.com/cms/cut\\_resistant\\_gloves](https://us.vwr.com/cms/cut_resistant_gloves).